

SOILTEST EQUIPMENTS



Auger Powerhead	85
Auger Hand Operated	85
Pocket Dial Penetrometer	86
Pocket Vane Shear Apparatus Set	86
Dynamic Cone Penetrometer	87
Speedy Moisture Tester	87
Universal Sample Extruder	88
Soil Lathe Device / Trimmer and Extruder	88
Laboratory Mixer	89
Mechanical End-Over-End Shaker	89
Semi Automatic Penetrometer for Liquid Limit	90
Digital Penetrometer for Liquid Limit	90
Manual Liquid Limit Device	91
Motorized Liquid Limit Device	91
Shrinkage Limit Test Set	92
Plastic Limit Test Set	92

Linear Shrinkage Mould	93
Munsell Soil Colour Chart	93
Front Loading Oedometer Test Set	94
Automatic Direct / Residual Shear Test Machine	95
Triaxial UU Test System Set	96
Automatic Soil Compactor	97
Standard Proctor Mould and Rammer	97
Modified Proctor Mould and Rammer	98
Automatic Cbr Test Machine	99-100
Cbr Mould and Accessories	101-102
In-situ CBR Test Apparatus Set	103
Sand Replacement Test Set	103
Sand Density Cone Test Set	104
Plate Bearing Test Set	104
Constant Head Permeability Set	105

www.testmak.com

AUGER POWERHEAD

Product Code

TMS-4800 | Auger Powerhead

TMS-4805 | Single Flute Bit Ø:80x1000 mm

TMS-4806 | Single Flute Bit Ø:100x1000 mm

TMS-4807 | Single Flute Bit Ø:150x1000 mm

TMS-4808 | Single Flute Bit Ø:200x1000 mm

Standards

ASTM D420, D1452 | AASHTO T86, T202 | CNR a. VI n.25



Desctription

Used in conjunction with sampling tubes to obtain disturbed or undisturbed soil samples. 47.7 cc displacement, 40:1 reduction ratio, 80-100-150-200 mm bit diameter. The Auger should be ordered separately.

P. Code	Dimensions(mm)	Weight (kg)	Power
TMS-4800	150X150X320 mm	12 kg	1600 W
TMS-4805	150x150x1150 mm	1 kg	1600 W
TMS-4806	200x150x1150 mm	2 kg	
TMS-4807	200x150x1150 mm	3 kg	1600 W
TMS-4808	200x150x1150 mm	4 kg	1600 W

HAND OPERATED AUGER SET

Product Code

TMS-4802 | Hand Operated Auger Set

TMS-4804 | Extension Rod for TMS-4802, 1 meter

Standards

ASTM D420, D1452 | AASHTO T86, T202 | CNR a. VI n.25



Desctription

Hand Auger and Sampling Set for augering down to a depth of 5 meters. The Hand Operated Auger Boring Set consists of a T-Handle with 1 m Rod, 80 mm, 100 mm and 150 mm diameter Auger Heads. Extension rods should be ordered separately

P. Code	Dimensions(mm)	Weight (kg)
TMC-3620	1100x200x200 mm	6 kg

POCKET PENETROMETER

Product Code

TMS-4810 | Pocket Dial Penetrometer, 0-6 kgf/cm²



Desctription

The Pocket Dial penetrometers are ideal instruments to determine the penetration resistance of cohesive soil, especially when various range measurements are required. Dual scale: 0-6 kgf/cm for 6.35 mm diameter plunger and 0-11 kgf/cm.

Technical Specifications	
Range (kgf/cm²)	0-6
Plunger Dia.	6.35 - 10 - 15 - 20 - 25
Dimensions(mm)	100x200x60 mm
Weight	0,5 kg

POCKET SHEAR VANE APPARATUS

www.testmak.com

Product Code

TMS-4325 | Pocket Vane Shear Apparatus Set

Standards

ASTM D 1558 | ASTM D 2573



Desctription

The Vane Shear Apparatus Set is used for determining resistance of soft clay, silt that are sensitive and cohesive soils at the place.

The Vane Shear Apparatus Set is supplied complete with:

- Standart Vane (0 -10 N /cm2) that is 25 mm
- Sensitive Vane Adaptor 0 -2 N /cm2
- High capacity vane adaptor 0 25 N/ cm2 and plastic mould.

Technical Specifications	
Dimensions(mm)	240x200x70 mm
Weight	0,5 kg

www.testmak.com

DYNAMIC CONE PENETROMETER

Product Code

TMS-4890 | Dynamic Cone Penetrometer (DCP)

Standards

BS 1377:9



Desctription

The Dynamic Cone Penetrometer is used for the rapid, in situ measurement of structural properties of existing road pavement constructed with unbound materials.

The design of the DCP is similar to that described incorporates an 8 kg weight dropping through a height of 575 mm and 60° cone having a diameter of 20 mm. with the standard DCP measurements can be made down to a depth of approximately 850 mm or when extension shafts are used to a recommended maximum depth of 2 m. Readings are usually taken after a set number of blows, changing the number according to the strength of the layer being penetrated. A typical test takes only a few minutes, therefore the instrument provides a very efficient method of obtaining information that would normally require the digging of test pits. The penetration hammer assembly consists of 8 kg hammer hammer shaft, anvil with plastic plate coupling for ruler and handle.

The Dynamic Cone Penetrometer is supplied complete with;

- · A hammer assembly,
- Penetration rod.
- 2 piece 60° cone.
- Metal plate coupling for ruler,
- · Segmented adaptor for extension rods,
- Segmented upper extension rod,
- · Segmented lower extension rod,
- 2 piece 13-17mm AF spanners.
- 3mm AF hex wrench,
- Tommy bar,
- A bottle of adhesive, 10cc.
- Users Manual
- Steel ruler.
- Carrying Case, heavy-duty wooden

SPEEDY MOISTURE TESTER

Product Code

TMS-4900 | Speedy Moisture Tester

Standards

BS 812-109



Desctription

Field conditions dust, pulp mixture and the soil, sand, clay, aggregate and other granular materials quickly and accurately detect moisture content. Moisture within the sample, as a result of reaction with calcium carbide acetylene gas occurs. The resulting pressure caused by gas, dial indicator is determined as the percent moisture content. Within a few minutes into the issue determines the moisture content of the sample. Precision Scale, measuring spoons and cups, calibration kit, 1 box of calcium carbide reagent, cleaning brush, wooden enclosure with a complete set as. 20 cup capacity Humidity Range: 0-20%.

The Speedy Moisture Tester is supplied complete with;

- · Vessel with Gauge
- · Digital Scale
- Scoop
- · Cleaning Brush
- Cleaning Cloth
- · Two Steel Pulverizing Balls
- · Plastic Case

Technical Specifications		
Capacity	20 Cup	
Humidity Range	0-20%	
Dimensions	510x380x200 mm	
Weight (kg)	8 kg	

UNIVERSAL SAMPLE EXTRUDER

Product Code

TMS-4350 | Universal Sample Extruder

Standards

EN 13286-2, 13286-47, 12697-30 | AASHTO T134, T180, T193, T245 | ASTM D698, D1557, D1883, D1559 | BS 1377-4, 1924-2, 598-107



Desctription

The TMS-4350 Sample Extruder is produced to easily extrude specimens from Marshall, CBR, standard and modified Proctor Moulds. The capacity of the extruder is 50 kN. Supplied complete with a manual hydraulic jack and 2 pcs. adaptor to extrude specimens from 100mm (4 "), 150 mm (6") inner diameter marshall, CBR standard and modified proctor, moulds.

TMS-4350 Technical Specifications		
Screw Travel 90 mm		
Ram Travel	130 mm	
Dimensions	300x300x500 mm	
Weight (kg)	30 kg	

SOIL LATHE DEVICE / TRIMMER AND EXTRUDER

www.testmak.com

Product Code

TMS-4860 | Soil Lathe Device / Trimmer and Extruder

TMS-4862 | Open Wire Saw TMS-4864 | Wire Saw TMS-4866 | Trimming Knife



Desctription

The Soil Lathe Device, Trimmer and Extruder is used to extrude and trim soil samples from 35 mm to 100 mm diameter to reduce samples. It should be used together with Wire Saw. Open Wire Saw, Wire Saw, Trimming Knife, Porcelain Mortar with Pestle.

Technical Specifications		
Specimen Lathe	35x70 mm to 100x200 mm	
Specimen Trimming and Extrusion	35x70 mm to 50x200 mm	
Dimensions	260 mm	

P. Code	Dimensions(mm)	Weight (kg)
TMS-4860	230x300x470 mm	16 kg

www.testmak.com

LABORATORY MIXER

Product Code

TMM-2430 | Laboratory Mixer - 220-240 V 50/60 Hz TMM-2430/110 | Laboratory Mixer - 110 V 60 Hz TMM-2421 | Stainles Steel Bowl 5 Liter Capacity

TMM-2422 | Whisk for TMM-2430

Standards

ASTM C109, C305 | AASHTO T106, T162 | EN 196-1, 196-3:2005, 413-2, 459-2, 480-1



Desctription

The Laboratory mixer is has been designed to mix mortars and cement pastes primarily to the requirements of standards. The mixing paddle revolves at a rate of 140 and 285 (r.p.m. clockwise) with a planetary motion of 62 and 125(r.p.m. clockwise) in low speed.

The Laboratory Mixer is supplied complete with;

- Stainles Steel Bowl, 5 L (approx.)
- Whisk for TMM-2430

Technical Specifications		
Dimensions	325x375x550 mm	
Weight (kg)	18 kg	
Power	400 W	

END-OVER-END SHAKER

Product Code

TMS-4870 | Mechanical End-Over-End Shaker, 220-240 V 50-60 Hz TMS-4870/110 | Mechanical End-Over-End Shaker, 110 V 60 Hz TMS-4872 | Gas Jar for TMS-4870

Standards

BS 1377:2



Desctription

The device is used for the determination of the particle density by the gas jar method and the particle size distribution by sedimentation. End-Over-End Shaker is used to rotate two 1 liter gas jars with rubber cover at 50 rpm. The 1 liter capacity gas jar is made of plexiglass and supplied complete with rubber cover. Gas Jar should be ordered separately.

P. Code	Dimensions(mm)	Weight (kg)	Power
TMS-4870	920x740x650 mm	24 kg	180 W
TMS-4872	120x120x270 mm	0,5 kg	

SEMI-AUTO PENETROMETER FOR LIQUID LIMIT

Product Code

Standards

TMS-4240 | Semi Automatic Penetrometer for Liquid Limit TMS-4240/110 | Semi Automatic Penetrometer for Liquid Limit TMS-4242 | 30° Penetration Cone

TMG-0568 | Sample Cups, Aluminium, Ø55 mm x h:35 mm

CEN ISO | TS 17892-12 | BS 1377:2 | NF P94-052-1





TMS-4242





DIGITAL PENETROMETER FOR LIQUID LIMIT

www.testmak.com

Product Code

TMS-4230 | Digital Penetrometer for Liquid Limit

TMS-4242 | 30° Penetration Cone

TMG-0568 | Sample Cups, Aluminium, Ø55 mm x h:35 mm

Standards

CEN ISO | TS 17892-12 | BS 1377:2 | NF P94-052-1



Desctription

The Semi-Automatic Liquid Limit Penetrometer test is based on the relationship between the moisture content at which clay soils pass from a plastic to a liquid state. The Soil Penetrometer is consists of steel base, leveling screws, digital penetration measurement gauge 0,01 mm precision, release button, automatic zeroing and spirit level.TMS-4240 is supplied complete with automatic timer unit. Penetration timer unit is used to release the plunger fitted with the needle to start the 5 seconds test.

Semi-Automatic Liquid Limit Penetrometer supplied complete with;

- · Automatic Penetration Timer Unit,
- 30° Penetration Cone
- Sample Cups, 3 pcs. Aluminium, Ø55 mm x h:35 mm.

P. Code	Dimensions(mm)	Weight (kg)
TMS-4240	220x300x410 mm	15 kg

Desctription

The Soil Penetrometer test is based on the relationship between the moisture content at which clay soils pass from a plastic to a liquid state. The Soil Penetrometer is consists of steel base, leveling screws, digital penetration measurement gauge 0,01 mm precision, release button, automatic zeroing and spirit level.

Digital Liquid Limit Penetrometer supplied complete with;

- 30° Penetration Cone
- Sample Cups, 3 pcs. Aluminium, Ø55 mm x h:35 mm.

P. Code	Dimensions(mm)	Weight (kg)
TMS-4230	220x170x410 mm	8

www.testmak.com

MANUAL LIQUID LIMIT DEVICE

Product Code

TMS-4250 | Manual Liquid Limit Device

TMS-4251 | Brass Cup

TMS-4252 | Brass Grooving Tool AASHTO

TMS-4253 | Metal Grooving Tool ASTM

TMS-4254 | Plastic Grooving Tool ASTM

TMS-4255 | Plastic Grooving Tool BS

Standards

Desctription

ASTM D4313 | AASHTO T89 | BS 1377:2







MOTORIZED LIQUID LIMIT DEVICE

Product Code

TMS-4260 | Motorized Liquid Limit Device

TMS-4251 | Brass Cup for TMS-4250 and TMS-4260

TMS-4252 | Brass Grooving Tool AASHTO

TMS-4253 | Metal Grooving Tool ASTM

TMS-4254 | Plastic Grooving Tool ASTM

TMS-4255 | Plastic Grooving Tool BS

Standards

ASTM D4313 | AASHTO T89 | BS 1377:2

Desctription

Motorized Liquid Limit Device is used to determine the moisture content Manual Liquid Limit Device is used to determine the moisture content at which clay soils pass from plastic to liquid state. The Liquid Limit Device consist of an adjustable crank and cam mechanism, a blow counter and a removable brass cup fitted on the base. The device supplied complete supplied complete with a brass grooving tool. with a brass grooving tool.

P. Code	Dimensions(mm)	Weight (kg)	
TMS-4250	200x220x150 mm	2 kg	

at which clay soils pass from plastic to liquid state. The Liquid Limit Device consist of an adjustable crank and cam mechanism, micro motor, a blow counter and a removable brass cup fitted on the base. The device

P. Code Dimensions(mm)		Weight (kg)	
TMS-4250	200x220x150 mm	2 kg	

SHRINKAGE LIMIT TEST SET

Product Code

TMS-4380 | Shrinkage Limit Test Set TMS-4381 | Shrinkage Prong Plate TMG-0104 | Porcelain Dish 120mm dia.

TMG-0569 | Aluminum Moisture Content Tin - Ø:45 mm h:10 mm TMG-0568 | Aluminum Moisture Content Tin - Ø:55 mm h:35 mm

TMG-0424 | Spatula 120 mm

TMG-0003 | Graduated Glass Cylinder 25 ml

Standards

ASTM D427 | ASTM T92



Desctription

Shrinkage Limit Test Set is used to determine the maximum moisture content at which the soil does not shrink after drying the sample.

The TMS-4380 Shrinkage Limit Test Set is supplied complete with:

- Porcelain Dish 120mm dia.
- Prong Plate
- Aluminum Moisture Content Tin Ø:45 mm h:10 mm, 2 pcs.
- Aluminum Moisture Content Tin Ø:55 mm h:35 mm
- · Spatula 120 mm
- Graduated Glass Cylinder 25 ml
- Carrying Case

P. Code	Dimensions(mm)	Weight (kg)
TMS-4380	350x300x100 mm	2 kg

PLASTIC LIMIT TEST SET

Product Code

TMS-4370 | Plastic Limit Test Set TMG-0992 | Glass Plate - 300x300x5 mm TMG-0104 | Porcelain Dish 120mm dia.

TMG-0570 | Aluminum Moisture Content Tin - Ø:75 mm h:30 mm

www.testmak.com

TMS-4372 | Steel Reference Rod Ø 3x100 mm

TMG-0424 | Spatula 120 mm

Standards

ASTM D4318 | AASHTO T90 | BS 1377:2



Desctription

Plastic Limit Test Set is determine as the lowest moisture content of a soil that will permit a sample to be rolled into threads of 3 mm diameter without the threads breaking.

The TMS-4370 Plastic Limit Test Set is supplied complete with:

- Glass Plate 300x300x5 mm
- Porcelain Mixing Dish Dia. 120 mm
- Moisture Content Tins Ø:75 mm h:30, 6 pcs.
- Steel Reference Rod
- · Spatula 120mm
- Carrying Case

P. Code	Dimensions(mm)	Weight (kg)
TMS-4370	360x360x150 mm	6 kg

www.testmak.com

LINEAR SHRINKAGE MOULD

Product Code

TMS-4460 | Linear Shrinkage Mould

Standards

BS 1377:2



Desctription

The Linear Shrinkage Mould is used to produce a specimen of 140 mm long x 12,5 mm radius. Shrinkage test covers the determination of linear shrinkage of soils and indicates the plastic properties of soils with a low clay content.

P. Code Dimensions(mm)		Weight (kg)	
TMS-4460	140x140x14 mm	0,4 kg	

HYDROMETER METHOD TEST KIT

Product Code

TMS-4300 | Hydrometer Method Test Kit

TMS-4301 | Heated Glass Water Tank - 600x300x380 mm

TMS-4302 | Mechanical Mixer 10.000 rpm - 220 V 50/60 Hz 1000W

TMS-4303 | Soil Hydrometer / 151H / 0.995 to 1038 g / ml 0.001 Precision

TMG-0008 | Graduated Glass Cylinders 1000 ml

TMG-0612 | Glass Thermometer 110 ° C


TMG-0015 | Beaker 250 ml



MUNSELL SOIL COLOR CHART

Product Code

TMS-4425 | Munsell Soil Colour Chart



Desctription

The Munsell Soil Chart provides a simple method for soil classification by of determining the color of soil specimens. Test set consists of 7 constant hue charts covering a total of 196 colors. The color chart and the diagram are fitted in a pocket size binder. Supplied complete with a Tropical Soil Color Chart.

P. Code Dimensions(mm)		Weight (kg)	
TMS-4425	155x200x35 mm	0,5 kg	

HYDROMETER METHOD TEST KIT

Standards

ASTM D422 | AASHTO T88 | UNE 103,102

Desctription

Hydrometer test kit is used to determine the quantitative such as silt and clay distribution of very fine particles in the ground. A heater and circulation motor allows a homogeneous of water temperature in water tank. Professional type mechanical mixer device is provided with set.

TMS-4300 Hydrometer Test Set consists of the following equipment;

- Heated Glass Water Tank 600x300x380 mm
- Mechanical mixer 10000 rpm
- Soil hydrometer / 151H / .995 to 1038 g / ml 0.001 Precision
- Graduated Glass Cylinders 1000 ml 6 pcs
- Glass Thermometer 110 ° C
- Beaker 250 ml
- P. Code
 Dimensions(mm)
 Weight (kg)

 TMS-4300
 400x900x350 mm
 15 kg

FRONT LOADING OEDOMETER TEST SET

Product Code

TMS-4110 | Front Loading Oedometer Test Set

TMS-4111 | Bench for Consolidation, 3 Oedometer Capacity

TMS-4112 | Front Loading Oedometer (Consolidation)

TMS-4113 | Consolidation Cell for High Pressure, Ø 50 mm

TMS-4114 | Consolidation Cell for High Pressure, ASTM Ø 63.5 mm (2.5")

TMS-4115 | Consolidation Cell for High Pressure, BS/EN, Ø 75 mm

TMS-4116 | Calibration Disc for Ø 50 mm Consolidation Cell

TMS-4117 | Calibration Disc for Ø 63,5 mm Consolidation Cell

TMS-4118 | Calibration disc for Ø 75 mm Consolidation Cell

TMS-4119 | Set of Weight for Consolidation - 80 kg (6x10 kg + 2x5 kg +

3x2 kg + 2x1 kg + 3x500 g + 2x250 g

TMS-4120 | Analog Dial Gauge, 30x0.01 mm

Standards

BS 1377 | ASTM D2435, D3877, D4546 | AASHTO T216

Desctription


Rigidly manufactured from aluminium alloy casting to provide a high degree of accuracy with any frame distorsion under load. The load bridge group is supported in high accuracy self-aligning seat balls. The beam provides three loading ratio: 9:1 10:1 11:1 and the beam assembly is fitted with an adjustable counterbalance weight. Maximun load: 150 kg of slotted weights. The oedometer accepts cells up to 100 cm2.

The TMS-4110 Front Loading Oedometer Test Set is supplied complete with:

- Bench for Consolidation, 3 Oedometer Capacity 1 pcs
- Front Loading Oedometer (Consolidation)
- Consolidation Cell for High Pressure, BS/EN, Ø 75 mm 3 pcs
- Calibration disc for Ø 75 mm consolidation cell 3 pcs
- Set of Weight for Consolidation 80 kg (6x10 kg + 2x5 kg + 3x2 kg + 2x1 kg + 3x500 g + 2x250 g) 3 set
- Analog Dial Gauge, 30x0.01 mm 3 pcs



www.testmak.com





www.testmak.com

DIRECT RESIDUAL SHEAR TEST MACHINE

Product Code

TMS-4180 | Automatic Direct / Residual Shear Test Machine, 220 V 50-60 Hz

TMS-4182-1 | Shear Box Assembly,60x60 mm

TMS-4182-2 | Specimen Cutter,60x60 mm

TMS-4182-3 | Extrusion Dolly, 60x60 mm

TMS-4183-1 | Shear Box Assembly,Ø 60 mm

TMS-4183-2 | Specimen Cutter,Ø 60 mm

TMS-4183-3 | Extrusion Dolly Ø 60 mm

TMS-4184-1 | Shear Box Assembly, 100x100 mm

TMS-4184-2 | Specimen Cutter, 100x100 mm

TMS-4184-3 | Extrusion Dolly 100x100 mm


TMS-4185-1 | Shear Box Assembly,Ø 100 mm

TMS-4185-2 | Specimen Cutter,Ø 100 mm

TMS-4185-3 | Extrusion Dolly Ø 100 mm

TMS-4186-1 | Shear Box Assembly, Ø 2.5 inch

TMS-4186-2 | Specimen Cutter,Ø 2.5 inch TMS-4186-3 | Extrusion Dolly Ø 2.5 inch





Standards

TMS-4181-6 | Slotted Weight Set,50 kg (4x10 kg + 1x5 kg + 2x2 kg + 1x1 kg) TS 1900-2 | ASTM D3080 | BS 1377-7 | AASHTO T236 | NF P094 071-1/2

Desctription

The System is capable of performing the consolidation and shearing phases of a standard direct shear and residual shear test under full automatic control. The system consists of a computer-controlled unit that utilizes micro-stepper motors to apply the vertical and horizontal loads to the soil specimen. The system is capable of doing a consolidation process at 32 phase automatically. Horizontal shearing can be applied at a specified rate of deformation or at a specified rate of horizontal force change, or at a specified set of force steps of

The system is capable of displaying the current status of a test and graphically portraying the progress of the test in real time. The system includes the capability for the operator to alter the test process and conditions at any stage of

Shear Box Assembly, Slotted Weight Set and other optional accessories including specimens cutter and extrusion be ordered separately.









TRIAXIAL TEST SYSTEM

Product Code

TMS-4700 | Triaxial UU Test System Set

TMS-4630 | Multi Speed Universal Electromechanic Test Machine

TMS-4624 | Triaxial Cell for 38 and 50 mm Samples

TMS-4625 | Triaxial Cell for 70 and 100 mm Samples

TMS-4710 | Oil and Water Constant Pressure System

TMS-4716 | De-Airring Water Tank with Hose

TMS-4765 | Block for Pressure Measurement

TMG-0916 | Pressure Transducer - 2000 kPa

Standards

ASTM D2850, D4767, D7181 | AASHTO T-297 | BS 1377-7, BS 1377-8

Desctription

The Triaxial UU Test System determining the mechanical properties of soils is a very important step to design foundations, embankments and other soil structures. The Triaxial UU Test System provides automated triaxial compression tests on cylindrical undisturbed and remolded soil samples. The stress-strain relation of undisturbed soil specimen can be investigated by subjecting the soil sample to different stress levels and drainage conditions with Triaxial Test System. De-Airring Water Tank, Block for Pressure Measurement, Pressure Transducer, Oil and Water Constant Pressure System should be ordered separately.


The TMS-4700 Triaxial UU Test System is supplied complete with:

- TMS-4630 | Multi Speed Universal Electromechanic Test Machine
- TMS-4624 | Triaxial Cell for 38 and 50 mm Samples
- TMS-4625 | Triaxial Cell for 70 and 100 mm Samples
- TMS-4701 | Software for Triaxial UU Test System

Technical Specification	
Maximium Load Capacity	50 kN
Infinitesimal Testing Speed	0,00001 - 6 mm/min
Speed Accuracy	± 0,5%
Vertical Daylight	0 - 790 mm
Distance Between Columns	360 mm
Base Diameter	158 mm
Software for Triaxial	Software for Triaxial UU Test System
Dimensions	430x380x1180 mm
Weight (approx)	105 kg



www.testmak.com



www.testmak.com

AUTOMATIC SOIL COMPACTOR

Product Code

TMS-4165 | Automatic Soil Compactor - 220-240 V 50-60 Hz TMS-4165/110 | Automatic Soil Compactor - 110 V 60 Hz

TMS-4166 | Compaction Rammer Face Ø 50.8 mm ASTM Standard TMS-4167 | Compaction Rammer Face Ø 50 mm EN/BS Standard

Standards

ASTM D558, D559, D560, D698, D1557, D1883 | EN 13286 2, 13286-47 | BS 1377:4 AASHTO T99, T134, T135, T136, T180, T193; NLT 107/98, 108/91, 111/87



Desctription

The Automatic Soil Compactor is designed to perform fast and accurate compaction of soil samples automatically to by international standards. The machine is designed to allow the hammer to drop the required height into the soil in the mould which rotates circularly to distribute the blows uniformly over the surface of the specimen in the mould. The drop height is adjustable to for EN Standard 300 mm and 450 mm, for ASTM Standard 305 mm and 457 mm. The compaction rammer is circular faced and interchangeable to for EN Standard 50 mm or for ASTM Standard 50.8 mm diameter. Compaction rammer weight is adjustable to 2.5 kg or 4.5 kg according to reference standard. The number of blows per layer can be set at the beginning of the compaction process by the digital counter according to the standard by user.

P. Code	Dimensions(mm)	Weight (kg)	Power
TMS-4165	430x250x1500 mm	140 kg	370 W

STANDARD PROCTOR MOULD AND RAMMER

Product Code

TMS-4445 | Standard Proctor Mould ASTM TMS-4446 | A Type Standard Proctor Mould EN TMS-4447 | Standard Proctor Mould BS

TMS-4450 | Standard Proctor Rammer ASTM TMS-4449 | A Type Standard Proctor Rammer EN TMS-4448 | 2,5 kg Standard Proctor Rammer BS

Standards

ASTM D 698, D 1557, D 558 | AASHTO T 99, T180, T 134 | EN 12386-2, 13286-4 | BS 1377:4,1924:2



Desctription

Standard Proctor Moulds are used for determining the moisture content and density of compacted soil. These are made of plated steel for corossion , includes collar, mould body and base plate. Standard Proctor Rammers are used to compact the soil sample in the Proctor. Different models are available conforming to the international standards.

P. Code	Internal Dia.(mm)	Body (mm)	Volume (cm³)	Weight (kg)
TMS-4445	101.6 ± 0.4 mm	116.4 ± 0.5 mm	944.0 ± 14	6,8 kg
TMS-4446	101.6 ± 0.4 mm	120 ± 1 mm	942	5 kg
TMS-4447	101.6 ± 0.4 mm	115,5± 0,5 mm	1	5,2 kg

P. Code	Internal Dia.(mm)	Body (mm)	Volume (g)	Weight (kg)
TMS-4450	50.8 mm	304.8± 1 mm	2495 ± 23 g	4.5 kg
TMS-4449	50 ± 0.5 mm	305± 3 mm	2500 ± 20 g	4.5 kg
TMS-4448	50 ± 0.5 mm	300± 3	2500 ± 25 g	4.5 kg

MODIFIED PROCTOR MOULD AND RAMMER

Product Code

TMS-4435 | Modified Proctor Mould ASTM TMS-4436 | B Type Modified Proctor Mould EN TMS-4437 | Modified Proctor Mould BS

TMS-4440 | Modified Proctor Rammer ASTM TMS-4439 | B Type Modified Rammer Mould EN TMS-4438 | 4,5 kg Modified Proctor Rammer BS



www.testmak.com

Standards

ASTM D 698, D 1557, D 558 | AASHTO T 99, T180, T 134 | EN 12386-2, 13286-4 | BS 1377:4,1924:2

Desctription

Modified Moulds are used for determining the moisture content and density of compacted soil. These are made of plated steel for corossion, includes collar, mould body and base plate. Modified Proctor Rammers are used to compact the soil sample in the Proctor. Different models are available conforming to the international standards.

P. Code	Desctription	Internal Dia.(mm)	Body (mm)	Volume (cm³)	Weight (kg)
TMS-4435	Modified Proctor Mould ASTM	152.4 ± 0.7 mm	116.4 ± 0.5 mm	2124 ± 25	10 kg
TMS-4436	B Type Modified Proctor Mould EN	150 ± 1 mm	120 ± 1 mm	2120	9,8 kg
TMS-4437	Modified Proctor Mould BS	152 ± 0.5 mm	127 ± 1 mm	2303	7,5 kg

P. Code	Desctription	Internal Dia.(mm)	Body (mm)	Volume (cm³)	Weight (kg)
TMS-4440	Modified Proctor Compaction Rammer ASTM	50.8 mm	457 ± 1.3 mm	4540 ± 10	8 kg
TMS-4439	B Type Modified Proctor Compaction Rammer EN	50 ± 0.5 mm	457 ± 3 mm	4500 ± 40	8 kg
TMS-4438	4.5 kg Compaction Rammer BS	50 ± 0.5 mm	450 ± 4 mm	4500 ± 50	8 kg

www.testmak.com

AUTOMATIC CBR TEST MACHINE

Product Code

TMS-4601 | Automatic Cbr Test Machine 50 kN - 220-240 V 50-60 Hz TMS-4601/110 | Automatic Cbr Test Machine 50 kN - 110 V 60 Hz

TMS-4164 | Linear Potentiometric Displacement Transducer, 25x0.001 mm

TMS-4470 | CBR Mould 6"

Standards

ASTM D1883 | BS 1377:4

Technical Specification		
Capacity	50 kN	
Platen Speed	0.625 to 5 mm/min	
Dimensions	400x560x1100 mm	
Weight (approx)	95 kg	
Power	370 W	



Desctription

The Automatic Cbr Test Machine is used to determine the maximum load and flow values of bituminous mixtures. Capacity is 50 kN. The machine comprises compact two column frame with adjustable upper cross beam. Platen speed is 0.625 to 5 mm/min. For safety, the up and down travel of the lower platen is limited the use of limit switches. Rapid adjustment of the platen is controlled using the up and down buttons on the front panel of the machine. The measuring system consists of a 50 kN capacity strain gauge load cell fitted to the upper cross beam to read stability values and the 25 x 0.001 mm linear potentiometric displacement transducer fitted to the breaking head. The Cbr machine is supplied complete with Cbr compaction mould.

Data acquisition and processing: by TCM readout unit featuring;

- Large graphic touch screen display 240 X 128 pixel,
- Effective resolution 16 bit,
- Effective sampling rate 40 Hz.
- Communication port

The Automatic Cbr Test Machine is supplied complete with;

- TMS-4601 | Automatic Cbr Load Frame
- TMS-4602 | Penetration Piston
- TMG-0906 | Load Cell, 50 kN
- TMG-0926 | Linear Potentiometric Displacement Transducer with Bracket, 25 x 0.001 mm
- TMS-4603 | Data acquisition and processing: by TCM LCD Display Unit Featuring
- TMS-4604 | PC Software
- TMS-4605 | Connection Cable



TCM304 LCD Graphic Display Data Acquisition and Control Unit

TCM304 LCD Graphic Display Data Acquisition and Control Unit is produced to control the machine and processing of data from load-cells, pressure transducers or displacement transducers which are fitted to the machine.

TCM304 LCD Graphic Display are controlled from the front panel consisting of a 240x120 pixel LCD display and function keys. One analog channel for load cell and one analog channel for displacement transducer exists.

TCM304 LCD Graphic Display are controlled with function keys on the front panel. One analog channel for load cell and one analog channel for displacement transducer exists.

MAIN FEATURES

- Automatically calculates flow and stability values
- 240*120 pixel blue-white graphic LCD display.
- High resolution 65.000 points.
- · Backlight function.
- · Capability of contrast calibration by light.
- 21 key touch membrane keyboard.
- Two analogical and two digital channels, use for load cell or pressure transducer etc.
- Standalone full automatic testing capacity.
- · Can make manual tests if requested.
- A sample type and dimension can be entered respect to the standards,
- Load-Time, Tensile-Time, Test Results and Sample reports observable and printable.
- One RS232 serial port for connecting either PC or printer for data transmission.
- · Comes with Connection cable and software.
- Large permanent memory up to 256 test results.
- Language select, English Turkish.

Data Acquisition & PC Software

CBR Test Software is developed for both EN 13286-47, BS 1377:4, ASTM D1883, AASHTO T193, NF P94-078 and UNI CNR 10009 standards CBR Tests. The Cbr Software includes control of machine, acquisition of load and displacement data, saving them and reports. The Software accepts specimen diameter and height as an input parameter. It automatically calculates correction factor coming from the standarts respect to specimen size. The stability value is calculated regarding to this factor. The software continously updates load and displacement until the end of test. When the test is completed, the sharpest slope of the graph is calculated. The point that this line crosses displacement axis is commented as an offset. This offset is subtracted from the displacement value at peak point and called as flow. Graphical outputs and reports can be saved as a MS Excel worksheet







100





www.testmak.com

www.testmak.com

CBR MOULD & ACCESSORIES - ASTM

Standards

ASTM D1883 | AASHTO T193

Desctription

Cbr Mould and Accessories is used for the laboratory evaluation of subgrade and subbase coarse materials in road construction.



P. Code	Desctription	Dimensions (mm)	Weight (kg)
TMS-4640	CBR Mould ASTM / AASHTO, plated steel mould body with 6" (152.4 mm) dia. x 7" (177.8 mm) height, supplied complete with extension collar and perforated base plate	200x200x450 mm	8,5 kg
TMS-4641	CBR Solid Base Plate ASTM - Plated Steel	200x200x100 mm	2 kg
TMS-4642	Perforated Base Plate - Plated Steel	200x200x100 mm	2 kg
TMS-4643	Spacer Disc with T handle ASTM, 150.8 mm dia. x 61.4 mm height"	200x200x300 mm	7 kg
TMS-4644	Annular Surcharge Weight ASTM, 2.27 kg	200x200x100 mm	2.3 kg
TMS-4645	Slotted Surcharge Weight ASTM, 2.27 kg	200x200x100 mm	2.3 kg
TMS-4646	Straightedge 300x30x5 mm	300x30x5 mm	1 kg
TMS-4647	Filter Paper for CBR Test No:5 x 150 mm dia. ASTM (Pack of 100)	200x200x20 mm	0,2 kg
TMS-4648	Filter screen, 144 mm dia. 150 µm mesh ASTM	150x150x20 mm	1 kg

CBR MOULD & ACCESSORIES - BS

Standards

BS 1377:4, 1924:2 | EN 13286-4

P. Code	Desctription	Dimensions (mm)	Weight (kg)
TMS-4650	CBR Type Mould BS (modified proctor) / Vibrating Hammer Mould BS, EN, TS-1900-1 internal dia.: 152 mm, internal height: 127 mm, supplied complete with extension collar and solid base plate	200x200x450 mm	8,5 kg
TMS-4651	Compaction plug with T handle BS, 150 mm dia. x 50 mm height	200x200x300 mm	7,5 kg
TMS-4652	Annular Surcharge Weight BS/EN, 2 kg	200x200x100 mm	2,3 kg
TMS-4653	Split Surcharge Weight BS/EN, 2 kg	200x200x100 mm	2,3 kg
TMS-4654	CBR Perforated Base Plate BS	200x200x100 mm	2 kg
TMS-4655	Filter Paper for CBR Test No:1 x 150 mm dia. BS (Pack of 100)	200x200x20 mm	0,2 kg
TMS-4656	C- Spanner for TMS-4650, 2 pcs	200x300x100 mm	1 kg
TMS-4657	Assembly tool for Base Plate for TMS-4650	350x30x15 mm	1 kg
TMS-4646	Straightedge 300x30x5 mm	300x30x5 mm	1 kg

CBR MOULD & ACCESSORIES - EN

Standards

EN 13286-4

Desctription

Cbr Mould and Accessories is used for the laboratory evaluation of subgrade and subbase coarse materials in road construction.



www.testmak.com

P. Code	Desctription	Dimensions (mm)	Weight (kg)
TMS-4660	B type Proctor Compaction Mould EN (modified), internal dia.: 150±1 mm, internal height: 120±1 mm, supplied complete with extension collar and solid base plate	200x200x180 mm	9 kg
TMS-4652	Annular Surcharge Weight BS/EN, 2 kg	200x200x100 mm	2,3 kg
TMS-4653	Split Surcharge Weight BS/EN, 2 kg	200x200x100 mm	2,3 kg
TMS-4661	CBR Perforated Base Plate EN	200x200x100 mm	2 kg
TMS-4662	Filter papers for CBR test, coarse, 148 mm dia.	200x200x20 mm	0.2 kg
TMS-4646	Straightedge 300x30x5 mm	300x30x5 mm	1 kg

SWEEL TEST EQUIPMENTS

Standards

EN 13286-47 | ASTM D1183 | AASHTO T193 | BS 1377:4, 1924:2

Desctription

Swell test equipment is used to monitor the swelling by placing it on top of the soil sample.



P. Code	Desctription	Dimensions (mm)	Weight (kg)
TMS-4671	Adjustable Stem and Perforated Plate for CBR Swelling Test	180x180x160 mm	2,5 kg
TMS-4672	Tripod for CBR Swelling Test	200x200x300 mm	1 kg
TMS-4673	Analog Gauge 30 mm travel x 0.01 mm division	150x100x80 mm	0,5 kg
TMS-4674	CBR Soaking Tank (6 pcs. CBR Mould Capacity)	600x600x400 mm	6 kg

www.testmak.com

IN-SITU CBR TEST APPARATUS SET

Product Code

TMS-4170 | In-situ CBR Test Apparatus Set

TMS-4171 | In-situ Test Frame

TMS-4172 | Load Ring 50 kN capacity

TMS-4173 | Penetration piston

TMS-4174 | Dial Gauge 30 mm travel x 0.01 mm with Holder

TMS-4175 | Set of extension rods (2 pcs. 102 mm, 1 pcs. 305 mm and 1 pcs. 610 mm length)

TMS-4176 | Datum bar assembly with two stands

Standards

ASTM D883 | D3668 | BS 1377:4, 1924 | AASHTO T193 | CNR UNI 10009 | NF P94-078



Desctription

In-situ CBR Test Apparatus is used to determine quickly and efficiently the bearing capacity of soils on road constructions, foundations, road subgrades. Load is applied through a mechanical jack and handwheel. Upper beam can be adjusted in height.

The In-situ CBR Test Apparatus Set is supplied complete with:

- Field Cbr Test Frame
- Load Ring 50 kN capacity
- Penetration piston
- Dial Gauge 30 mm travel x 0.01 mm with Holder
- \bullet Set of extension rods (2 pcs. 102 mm, 1 pcs. 305 mm and 1 pcs. 610 mm length)
- Datum bar assembly with two stands
- Sloted Surcharge Weight 2.27 kg, 2 pcs
- Anular Surcharge Weight 4.54 kg
- Anular Surcharge Weight 9.08 kg
- · Vehicle bracket and wooden carrying case

IN-SITU CBR TEST APPARATUS SET

Product Code

TMS-4130 | Sand Replacement Test Set 100 mm

TMS-4131 | Sand Pouring Cylinder for TMS-4130

TMS-4132 | Calibration Container for TMS-4130

TMS-4133 | Metal Tray for TMS-4130

TMS-4134 | Sand Replacement Test Set 150 mm

TMS-4135 | Sand Pouring Cylinder for TMS-4134

TMS-4136 | Calibration Container for TMS-4134

TMS-4137 | Metal Tray for TMS-4134

TMS-4138 | Sand Replacement Test Set 200 mm

TMS-4139 | Sand Pouring Cylinder for TMS-4138

TMS-4140 | Calibration Container for TMS-4138

TMS-4142 | Metal Tray for TMS-4138

Standards

BS 1377:9, 1924:2



Desctription

Sand Replacement Test Set is used to determine the in-situ density of fine grained compacted soil. The test consists in digging a hole into the ground and then collect, dry and weight the sampled soil. The hole is than filled with dry sand from the cone container.

P. Code	Dimensions(mm)	Weight (kg)
TMC-3430	300x300x450 mm	8 kg
TMC-3434	300x300x500 mm	15 kg
TMC-3438	500x500x660 mm	26 kg

IN-SITU CBR TEST APPARATUS SET

Product Code

TMS-4145 | Sand Density Cone Test Set 6,5"

TMS-4146 | Double Cone

TMS-4147 | Metaş Tray for TMS-4145

TMS-4148 | Plastic Sand Jar 5 liter

Standards

ASTM D1556 | AASHTO T181, T191



Desctription

Sand Density Cone Test Set is used to determine the in-situ density of fine grained compacted soil. Includes double cone, plastic sand jar and metal tray with flanged hole.

P. Code	Dimensions(mm)	Weight (kg)
TMC-4145	300x300x550 mm	3,5 kg

PLATE LOADING TEST SET

Product Code

TMS-4220 | Plate Bearing Test Set 200 kN Capacity TMS-4222 | Plate Bearing Test Set 500 kN Capacity

Standards

ASTM D 1194 | ASTM D 1195 | ASTM D 1196 | BS 1377



www.testmak.com

Desctription

Plate Bearing Test Set is performed for the determination of the bearing capacity of a soil in-situ on road constructions, foundations, road subgrades, airport and highway pavements. A wide range of plate bearing test equipment are available, together with many accessories according to the different Standards and specific enduser needs. All test sets supplied complete with 1,5 m long flexiable hose with quick release coupling.

The TMS-4220 Plate Bearing Test Set 200 kN is supplied complete with:

- 200 kN Capacity Piston Assemly
- Hydraulic Hand Pump,
- Pressure Transducer,
- Datum Bar 2.4 m long , 3 pcs.
- 25 mm travel x 0.01 mm digital dial gauges with dial supports,
- 300 mm and 450 mm dia. loading plates,

The TMS-4222 Plate Bearing Test Set 500 kN is supplied complete with:

- 500 kN Capacity Piston Assemly
- Hydraulic Hand Pump,
- Pressure Transducer,
- Datum Bar 2.4 m long , 3 pcs.
- 25 mm travel x 0.01 mm digital dial gauges with dial supports,
- 600 mm and 760 mm dia. loading plates

P. Code	Dimensions(mm)	Weight (kg)
TMC-4220	650x350x600 mm	130 kg
TMC-4222	850x850x130 mm	150 kg

www.testmak.com

CONSTANT HEAD PERMEABILITY SET

Product Code

TMS-4390 | Constant Head Permeability Set for Ø 75 mm cell

TMS-4391 | Constant Head Permeability Set for Ø 114 mm cell

TMS-4392 | Constant Head Permeability Cell, 75 mm dia.

TMS-4393 | Constant Head Permeability Cell, 114 mm dia.

TMS-4394 | Wooden Stand with 3 Manometer Tubes

TMS-4395 | Constant Level Water Tank, 7 Liter



BS 1377 | ASTM 46-47



Desctription

Constant Head Permeability Sets are used to determine the permeability of granular, gravel and sand soils. The specimen is formed in an acrylic permeability cell, and water is passed through it from a constant level tank. The permeability cell has pressure points at different levels which are connected to the manometer tubes fixed on a stand with graduated scale. Two constant head permeability cells are available: 75 mm and 114 mm diameter. Constant Level Water Tank Should be ordered separately.

The TMS-4390 Constant Head Permeability Set is supplied complete with;

- TMS-4392 | Constant Head Permeability Cell, 75 mm dia.
- TMS-4394 | Wooden Stand 1500 mm with 3 Manometer Tubes 1000 mm long
- TMS-4396 | 3 m Hose

The TMS-4391 Constant Head Permeability Set is supplied complete with;

- TMS-4393 | Constant Head Permeability Cell, 114 mm dia.
- TMS-4394 | Wooden Stand 1500 mm with 3 Manometer Tubes 1000 mm long
- TMS-4396 | 3 m Hose

P. Code	Dimensions(mm)	Weight (kg)
TMC-4390	220x70x1500 mm	9,5 kg

P. Code	Dimensions(mm)	Weight (kg)
TMC-4391	220x70x1500 mm	13 kg

FALLING HEAD PERMEABILITY SET

Product Code

TMS-4401 | Falling Head Permeability Cell 100 mm dia.

TMS-4402 | Wooden Stand with 4 Manometer Tubes

TMS-4403 | Soaking Reservoir Tank

TMS-4404 | 3 m Hose

TMS-4405 | De-Airing tank with in -let and a flow outlet Connection

TMS-4406 | Single Stage Vacuum Pump

Standards

BS 1377 | ASTM 46-47

Desctription



www.testmak.com

The Falling Head Permeability Set is used to study the behaviour of soil, particularly fine grained soils such as clay-like or silty soils, with respect to water flow. The Falling Head Permeability Cell is produced from plated steel with an inside diameter of 100 mm. The Wooden Stand is fitted with 4 glass Manometer Tubes of each 1500 mm long with inside diameters of about 21 mm, 12 mm, 5 mm and 3.5 mm. All tubes have connection valves. The Soaking Reservoir Tank is produced from plated steel with an over-flow tube and is used for containing the permeability cell during the test. De-Airing tank with in -let and a flow outlet Connection and Single Stage Vacuum Pump Should be ordered separately.

The TMS-4400 Constant Head Permeability Set is supplied complete with;

- TMS-4401 | Falling Head Permeability Cell 100 mm dia.
- TMS-4402 | Wooden Stand with 4 Manometer Tubes
- TMS-4403 | Soaking Reservoir Tank
- TMS-4404 | 3 m Hose

P. Code	Dimensions(mm)	Weight (kg)
TMC-4390	240x100x1500 mm	14 kg